


Create Alert(s)

1. Click the "Alerts" icon on the toolbar
2. Click on "Add New Alert" button
3. Click on "Add Device" button
4. Type in the alert name
5. Select the alert type from the drop-down box
6. Fill out all required information. (Every alarm type has different fields)
7. Check the information and click on "Save" button



Create Alert

5 Alert Type: Account Wide: ☐ Yes ☒ No

3 Add Devices: 0 device(s) selected - view selected devices

4 Alert Name*:

Every*: Hours

Scheduled For: Except on:

6

Send Alert to:

Email

add another

7

